

MAX Series

C-arm Table by Surgical Tables, Inc



Easy c-arm access

500 lb weight capacity

2" patient comfort pad

Light weight and easy to move

24" wide top

Made in the USA

Standard Accessories:

Hand Control




Patient Straps



Optional Accessories



Standard Side Rails



Paper Roll



Foot Control

*Streamline 1, 2, & 3 require upgrade to use both hand and foot controls simultaneously



Modular Side Rails



Surgical Arm Boards

*Require side rails



Procedural Arm Boards

Configurations/Specifications

MAX-ERCP Series

MAX4-ERCP

Height 30"-40"
Lateral Tilt 17°
Trendelenberg 13°
Longitudinal Travel 16"


MAX5-ERCP

Height 33.5"-43.5"
Lateral Tilt 17°
Trendelenberg 13°
Longitudinal Travel 16"
Lateral Travel 10"

Standard MAX Series


MAX 1

Height 28"-38"




MAX 2

Height 28"-38"




MAX 3

Height 28"-38"




MAX 4

Height 30"-40"



MAX 5

Height 33.5"-43.5"



Surgical Tables, Inc

2 DeBush Ave Unit C3
Middleton, MA 01949 USA
888-737-5044/978-777-4031
sales@surgicaltables.com
www.surgicaltables.com

500 lb weight limit
Tabletop: 84"L x 24"W
120V/60Hz or 230V/50Hz
5" Locking Swivel Castors
2" Thick Mattress
Lateral Tilt: 17°
Trendelenberg: 13°
Longitudinal Travel 10"
Lateral Travel: 10"
Table Dimensions: 86"L x 26"W