


VMAX-Pro Series



Vascular 4-way float C-arm Table
by Surgical Tables, Inc



Instant-Response 4-way float

Easy c-arm access

500 lb weight capacity

2" patient comfort pad

Light weight and easy to move

24" wide top

Made in the USA

Standard Accessories:

Joystick



Patient Straps



Fixed Side Rails



Optional Accessories



Procedural Arm Boards



Paper Roll



Foot Control



Catheter Tray



Modular Side Rails



Surgical Arm Boards
*Require siderails



Rectangular Arm Board



Lollipop Arm Board



*Trendelenberg only available for VMAXT-Pro

Configurations/Specifications



VMAX-Pro

Height 28"-38"
4 Way Float
Longitudinal Travel 36" (+/- 18")
Lateral Travel 14" (+/- 7")

VMAXT-Pro

Height 30"-40"
4 Way Float
Longitudinal Travel 36" (+/- 18")
Lateral Travel 14" (+/- 7")
Trendelenberg +/- 15°



Surgical Tables, Inc

2 DeBush Ave Unit C3
Middleton, MA 01949 USA
888-737-5044/978-777-4031
sales@surgicaltables.com
www.surgicaltables.com

500 lb weight limit
Tabletop: 84"L x 24"W
120V/60Hz
5" Locking Swivel Castors
2" Thick Mattress
Trendelenberg: +/- 15°
Longitudinal Travel 36" (+/- 18")
Lateral Travel: 14" (+/- 7")
Table Dimensions: 86"L x 26"W